Elaboro GmbH
Hagenower Str. 73
19061 Schwerin | Germany

TEL : +49 385 3993 300
E-MAIL : info@elaboro.de

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CERAFUSION® LiSi SCAN | MATTIFY - Scan spray for ultra fine details

1. Application Area

CERAFUSION® LiSi SCAN is a high-detail coating spray for the temporary matting of surfaces used in optical measuring and scanning procedures. It provides reliable matting even on highly demanding surfaces and, thanks to a particle size of only 1-2 µm, enables extremely precise detail reproduction.

The product can be applied to a wide range of materials, including plaster, metals, plastics, glass, and ceramics. It also delivers a uniform and fine coating on greasy or slightly moist surfaces. After application, the layer can be removed without residue using water or steam. The spray is free of titanium dioxide, highly economical in use, and odorless. Intraoral use is not permitted.

2. General Safety & Handling Instructions

CERAFUSION® LiSi SCAN is intended exclusively for use by trained professionals. Please observe the instructions on the spray can:

- Flammable aerosol. Keep out of reach of children.
- Pressurized container - may burst if heated.
- Keep away from heat, hot surfaces, sparks, and open flames. Do not smoke.
- Protect from sunlight. Do not expose to temperatures exceeding 50 °C / 122 °F.
- Do not pierce or burn, even after use.
- Use only in well-ventilated areas. Do not inhale spray mist.
- Wear protective clothing, eye protection, and a mask; use workplace extraction if necessary.
- Causes serious eye irritation. Seek medical advice if irritation persists.
- Do not force open or incinerate after use.
- Always empty spray cans completely.
- Intraoral use is not permitted. Follow the safety data sheet.

3. Storage and Working Conditions

- Store at room temperature, max. -20 °C to 45 °C / -4 °F to 113 °F
- Processing temperature: 15-35 °C / 59-95 °F
- Excessively high or low temperatures impair spray performance
- Use only in well-ventilated rooms
- Use extraction systems and dust protection masks
- Ensure good lighting to identify uniform coating

4. Application

Preparation

- Shake the can vigorously before use until the mixing balls are audible.
- Attach the spray head with nozzle and shake again briefly.
- Shake again after each spraying pause to prevent nozzle clogging.

Test & Spraying Technique

- Spray distance: approx. 15-20 cm from the object.
- Apply in short spray bursts.
- Hold the can as upright as possible.
- Use only the supplied spray head.

Notes for Optimal Application

- Wet spots indicate that the spray distance is too short.
- Uneven powder deposits indicate overdosing or insufficient shaking.
- If the spray head is clogged, clean it with compressed air. Alternatively, it can be cleaned in an ultrasonic bath.

Explanation of danger and warning signs

Manufacturer

Follow our
instructions for use!


Caution,
observe warning notices


Item No.

Charge No.

Working
Temperature

Transport & Storage
Temperature

Expiry Date

Caution! Flammable Aerosol
Causes severe eye irritation!

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©2019-2026, Elaboro GmbH