Elaboro GmbH
Hagenower Str. 73
19061 Schwerin | Germany

TEL : +49 385 3993 300
E-MAIL : info@elaboro.de

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CERAFUSION® LiSi CONDITIONER | CONNECT - Pre-Bond spray for zirconia

Raw ceramic for the manufacture of dental prostheses. Type, classification, and labeling according to DIN EN ISO 6872: Dental ceramics type I, class 1

1. Area of application

LiSi CONDITIONER is a ready-to-use, easy-to-apply lithium silicate glass ceramic for the adhesive bonding of monolithic zirconia restorations. The surface acquires glass ceramic properties, allowing crowns, bridges, veneers, inlays, onlays, Marylands, etc. to be bonded directly like glass ceramics.

No pre-treatment of the surface is necessary—no bonders, solvents, or pre-firing are required.

The zirconia surfaces should be free of dust and grease to ensure even coating.

LiSi CONDITIONER is suitable for all dental zirconia.

2. General safety and work instructions

LiSi CONDITIONER is intended exclusively for use in dental laboratories by trained specialists. Please follow the instructions on the spray can:

- Flammable aerosol, keep out of reach of children
- Container is pressurized - may burst if heated
- Keep away from heat, hot surfaces, sparks, and open flames, no smoking
- Protect from sunlight; do not expose to temperatures exceeding 50°C / 122°F
- Do not pierce or burn. (even after use)
- Use only in well-ventilated areas, do not inhale spray mist
- Wear protective gloves, protective clothing, eye protection, and a mask, use workplace extraction if necessary
- Causes serious eye irritation; if eye irritation persists, seek medical advice
- Do not force open or burn after use
- Always empty spray cans completely
- Intraoral use is not permitted, refer to the safety data sheet

3. Storage and working conditions

- Store at room temperature
- Processing temperature: 15-35 °C / 59-95 °F
- Temperatures that are too high or too low will impair the spray result
- Use only in well-ventilated areas
- Use an extraction system and dust mask
- Ensure good lighting to detect uneven coating

4. Preparation of the zirconia restoration

The finished sintered and fitted zirconia restoration must be dry, clean, and free of dust and grease in order to achieve optimal results. A very thin layer of LiSi CONDITIONER is sufficient.

For fully anatomical monolithic restorations, functional and aesthetic features such as chewing surfaces, contact points, color gradients, and effects should be taken into account before sintering the zirconia. Crown margins, marginal fits, and occlusion are practically unaffected by a single layer application.

We recommend applying the final layer of LiSi CONDITIONER only after dental try-in and any necessary corrections.

5. Application

Preparation

- Shake the spray can vigorously for at least 1 minute before use until the mixing balls can be heard
- Attach the spray head with nozzle and shake the can again briefly
- Shake again after each spray pause to prevent nozzle clogging

Testing & spraying technique

- Before using for the first time, test spray (e.g., on glass)
- Spray distance: approx. 15-20 cm from the object
- **Apply in short bursts**
- Hold the spray can as upright as possible
- Only use the spray nozzle provided

Tips for optimal use

- Wet spots indicate that the spray distance is too short
- Uneven powder deposits indicate overdosing or insufficient shaking
- Always apply LiSi CONDITIONER in short bursts
- Only spray enough so that the zirconium surface still shimmers through the white powder

6. Troubleshooting before ceramic firing

Incorrectly sprayed parts can be easily washed off with water or steamed off. The same applies if too much material has been applied. After drying with compressed air, LiSi CONDITIONER can be sprayed on again. If small areas are damaged when handling the sprayed dental restorations, these can be easily sprayed over again with a thin layer.

7. Cleaning

Clean the spray head immediately after using the spray can (e.g., with lukewarm water and compressed air, or an ultrasonic bath if necessary). Then dry with compressed air.

8. Instructions for firing

Perform the ceramic firing according to the specifications, see table diffusion firing. During the diffusion firing, components of the zirconium oxide and components of LiSi CONDITIONER diffuse into each other in the area close to the surface and form a strong bond during cooling. At the same time, a smooth surface is formed.

a. Diffusion firing: When applied correctly, you will obtain a resilient, smooth, and homogeneous surface of the highest quality after the first firing, with optimal adhesion to the zirconium oxide. The diffusion firing of LiSi CONDITIONER can be carried out simultaneously with LiSi PURE.

b. Alternative firing Recommendations: For older ceramic ovens, you can alternatively use the simplified firing program without vacuum.

9. Surface conditioning for adhesive bonding

a. Pretreatment of the coated restoration by etching

The classic bonding protocol for glass ceramics applies:

- Application of a ceramic etching gel (e.g., IPS Ceramic Etching Gel)
- Etching time: 2-3 minutes (follow manufacturer's instructions)

Etching creates adhesive surfaces on the all-ceramic restoration and strengthens the adhesive effect between the composite and the ceramic surface. In most cases, the conditioned surface still appears mirror-smooth to the naked eye despite etching.

The manufacturer's instructions for use must be strictly adhered to.

b. Adhesive Bonding

For the adhesive bonding of zirconia restorations coated with LiSi CONDITIONER, we recommend:

- Universal single-component bonding agents for creating a bond between the fixing composite and glass or oxide ceramics
e.g., Monobond Plus (application time: 60 seconds)

For permanent bonding, Variolink Esthetic from Ivoclar is suitable, among other products.

10. Tips for adjusting the vacuum setting

The vacuum used must not exceed 30%. The corresponding values are listed in the firing parameters table.

If the vacuum on the device cannot be freely adjusted, do not use vacuum at all.

11. Troubleshooting after ceramic firing

The LiSi CONDITIONER coating is only 5-10µm thick when applied normally. If small pores are visible on the surface after firing, this is due to insufficient powder application. Spray the restoration again—without any surface treatment—and repeat the firing process. If necessary, check the oven calibration and firing parameters.

DIFFUSION FIRING												Vacuum					
Closing Time	Standby Temperature		Heating Rate		Firing Temperature		Holding Time	Heating Rate		Firing Temperature		Holding Time	Long-term Cooling	Cooling Rate		Vacuum	
min	°C	°F	°C/min	°F/min	°C	°F	min	°C/min	°F/min	°C	°F	min	depends on size	°C/min	°F/min	400°C	752°F
1	400	752	40*	72*	820	1508	5	20	36	920	1688	5	yes	20 - 80	35 - 140	max. 30%	

ALTERNATIVE FIRING RECOMMENDATIONS FOR OLDER CERAMIC OVENS											
Closing Time	Standby Temperature		Heating Rate		Firing Temperature		Holding Time	Long-term Cooling	Cooling Rate		Vacuum
min	°C	°F	°C/min	°F/min	°C	°F	min	depends on size	°C/min	°F/min	-
1	400	752	40*	72*	920	1688	5	yes	20 - 80	35 - 140	no

*For veneers and individual crowns, a heating rate of 80°C/min · 162°F/min can be used.

Explanation of danger and warning signs

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©2015-2026, Elaboro GmbH