Elaboro GmbH
Hagenower Str. 73
19061 Schwerin | Germany

TEL : +49 385 3993 300
E-MAIL : info@elaboro.de

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CERAFUSION® LiSi LOWFUSE | REFINE - Glaze & repair spray for ceramics

Raw ceramic for the manufacture of dental prostheses. Type, classification, and labeling according to DIN EN ISO 6872: Dental ceramics type I, class 1

1. Area of application

LiSi LOWFUSE is the preferred choice as a final finish and sealing layer for zirconia restorations refined with LiSi PURE in combination with stains to achieve the highest aesthetic results in individual work.

LiSi LOWFUSE is a ready-to-use lithium silicate thin-film glaze for surface finishing of dental ceramics. It is used to finalize ceramic restorations such as crowns, bridges, or monolithic work to create a smooth, pore-free, and high-gloss surface. It is suitable for glass and pressable ceramics, as well as for sealing stains and optimizing already glazed surfaces without changing the fit of the restoration.

Due to the low firing range of around 800 °C, surfaces can be refined, repaired, or sealed, eliminating the need for polishing while improving plaque resistance and aesthetics.

2. General safety and work instructions

LiSi LOWFUSE is intended exclusively for use in dental laboratories by trained specialists. Please observe the instructions on the spray can:

- Flammable aerosol, keep out of reach of children
- Container is pressurized - may burst if heated
- Keep away from heat, hot surfaces, sparks, and open flames, no smoking
- Protect from sunlight; do not expose to temperatures exceeding 50°C / 122°F
- Do not pierce or burn. (even after use)
- Use only in well-ventilated areas, do not inhale spray mist
- Wear protective gloves, protective clothing, eye protection, and a mask, use workplace extraction if necessary
- Causes serious eye irritation; if eye irritation persists, seek medical advice.
- Do not force open or burn after use
- Always empty spray cans completely
- Intraoral use is not permitted, refer to the safety data sheet

3. Storage and working conditions

- Store at room temperature
- Processing temperature: 15-35 °C / 59-95 °F
- Temperatures that are too high or too low will impair the spray result
- Use only in well-ventilated areas
- Use an extraction system and dust mask
- Ensure good lighting to detect uneven coating

4. Preparation of the restoration

The finished, polished, and adjusted restoration must be dry, clean, and free of dust and grease in order to achieve optimal results. Before application, polishing paste, grinding, and other contamination residues must be completely removed, for example by suitable fine blasting and subsequent steam blasting or cleaning in an ultrasonic bath. Bare zirconia must be coated with LiSi PURE before application to create a suitable adhesive base for the subsequent glaze. Glazing with LiSi LOWFUSE requires only a small amount of material; the spray intensity must be adjusted according to the porosity and roughness of the ceramic. For fully anatomical monolithic restorations, functional and aesthetic features such as chewing surfaces, contact points, and color and effect designs should be taken into account before glazing. Crown margins, marginal fit, and occlusion are not affected by a single layer application (10-20µm).

5. Application

Preparation

- Shake the spray can vigorously for at least 1 minute before use until the mixing balls can be heard
- Attach the spray head with the lance and shake the can again briefly
- Shake again after each spraying interval to prevent the nozzle from clogging

Testing & spraying technique

- Before using for the first time, test spray (e.g., on glass)
- Spray distance: approx. 15-20 cm from the object
- **Apply in short bursts**
- Hold the spray can as upright as possible
- Only use the spray nozzle provided

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- Wet spots or streaks indicate that the spray distance is too short.
- Uneven powder deposits indicate overdosing or insufficient shaking
- Always apply LiSi LOWFUSE in short bursts
- The spray should only be applied so heavily that the restoration still shines through the white powder slightly

6. Troubleshooting before ceramic firing

Incorrectly sprayed parts can be easily washed off with water or steamed off; the same applies if too much material has been applied. After drying with compressed air, LiSi LOWFUSE can be sprayed on again. If small areas are damaged when handling the sprayed dental restorations, these can be easily sprayed over again with a thin coat.

7. Tips for firing

Perform the ceramic firing according to the specifications, see table Firing parameters. During the diffusion firing, components of the ceramic and components of LiSi LOWFUSE diffuse into each other in the area close to the surface and form a strong adhesive bond during cooling. At the same time, a smooth surface with a layer thickness of 10-20µm is formed.

Glaze firing: When applied optimally, you will obtain a resilient, smooth, and homogeneous surface of the highest quality after the first firing. Depending on the surface texture of the dental ceramic, it may be necessary to spray and fire another layer. Firing can be carried out in combination with high-melting stains and ceramic materials. For older ceramic ovens, you can alternatively use the firing program without vacuum.

8. Tips for adjusting the vacuum setting

The vacuum used must not exceed 30%. The corresponding values are listed in the combustion parameters table. If the vacuum on the device cannot be freely adjusted, vacuum must not be used at all.

9. Troubleshooting after ceramic firing

The LiSi LOWFUSE coating is 10-20µm thin under normal application. If small pores are visible on the surface after firing, this is due to insufficient powder application. Spray the restoration again—without any surface treatment—and repeat the firing process. If necessary, check the oven calibration and firing parameters.

10. Cleaning

Clean the spray head immediately after using the spray can (e.g., with lukewarm water and compressed air, or an ultrasonic bath if necessary). Then dry with compressed air.

Tips for the practitioner

During professional teeth cleaning (PTC), surfaces treated with LiSi LOWFUSE must not be blasted with powder materials. Roughening of the glass ceramic is irreversible and can increase susceptibility to plaque. Please also note the information on workplace exposure, transport, and storage in the safety data sheet.

GLAZE FIRING											
Closing Time	Standby Temperature		Heating Rate		Firing Temperature		Holding Time	Long-term Cooling	Cooling Rate		Vacuum
min	°C	°F	°C/min	°F/min	°C	°F	min	depends on size	°C/min	°F/min	-
1	400	752	50	90	800	1472	3	yes	20 - 80	35 - 140	no

Explanation of danger and warning signs

Follow our instructions for use!


Caution, observe warning notices


Item No.

Charge No.

Working Temperature

Transport & Storage Temperature

Expiry Date

Caution! Flammable Aerosol Causes severe eye irritation!

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©2015-2026, Elaboro GmbH